

New!

Distributed by:



217 Industrial Court, Wabasha, MN 55981

1-800-363-7846

www.medi-stim.com

STIM-U-WEAR™

Conductive Garment Electrodes

Stim-U-Wear conductive garments offer comfortable electrotherapy stimulation to areas that are not effectively treated by traditional adhesive electrodes. These form-fitting garments made of highly conductive antibacterial silver yarns, provide stimulation over the entire surface area of the garment, eliminating hot spots often found in other garments. This is an effective treatment option for pain management or muscle rehabilitation.

- ★ Effective treatment for patients suffering from Carpal Tunnel, Arthritic Pain, Diabetic Neuropathy, and Repetitive Strain Disorders
- ★ Increases circulation, reduces edema and post-operative swelling and pain
- ★ Provides comfortable, even stimulation throughout the entire surface area of the garment electrode
- ★ Compatible with all makes and models of electrotherapy stimulators (TENS, NMES, IF, HVPC, MENS)
- ★ Non-gelled alternative electrode for patients with skin sensitivities
- ★ Reimbursed by Medicare and most insurance companies
- ★ High-end quality, low end price!
- ★ HCPC Code: E0731



Knee Sleeve

Part #CSL-L

Large (5.5" W x 10.5" L)

Max. Circ. 25"



Elbow Sleeve

Part #CSL-S

Small (4.5" W x 7" L)

Max. Circ. 15.5"



Sock

Part #CSK-0

One Size Fits Most

(3.5 W x 13" L)

Distributed by:

Glove

Part #CGL-S

Small (3" W x 7" L)

Part #CGL-L

Large (3.5" W x 8.5" L)



Garments are sold individually, not in pairs, and do not include a return electrode that is necessary for the garment to conduct the current. All garments may have a +/- 2cm tolerance of sizes. Composition: 76% Nylon, 16% Silver yarn, 8% Spandex. Connector: 3.5mm snap + 2mm pin socket.

