



# BioStim<sup>®</sup> LX

## TENS UNIT

### Features

- ▶ Dual Channel (4-lead) All Digital Tens
- ▶ Four Treatment Modes
- ▶ Flip-top Panel Design
- ▶ Removable Plastic Belt Clip
- ▶ Power ON/OFF Button
- ▶ Treatment Timer
- ▶ 5-Year Warranty

### Technical Specifications

- ▶ Channel Dual
- ▶ Amplitude Continuously adjustable from 0 - 98 mA peak
- ▶ Waveform Asymmetrical, biphasic square wave
- ▶ Pulse Rate 1 - 150 Hz
- ▶ Pulse Width 10 - 250  $\mu$ S, adjustable
- ▶ Modes: Continuous, Burst 1, Burst 2, Modulation
  - *Continuous*: Continuous stimulation
  - *Burst 1*: 0.5 seconds on, 0.5 seconds off. Adjustable pulse rate from 1 - 150 Hz
  - *Burst 2*: 8 pulses per burst / 2 bursts per second
  - *Modulation*: Decrease / Increase of set width of 50% over a 5 second period
- ▶ Power Supply 4 AA batteries
- ▶ Size 8.2 x 7 x 4.5 cm
- ▶ Weight 226 grams



Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 48" pin lead wires, 4 each AA alkaline batteries, zippered case pouch, and instruction manual

Distributed by:

Medi-Stim, Inc., 217 Industrial Court  
Wabasha, MN 55981 USA  
Ph. 1-800-363-7846 Fax 1-651-565-2410  
[www.medi-stim.com](http://www.medi-stim.com) [sales@medi-stim.com](mailto:sales@medi-stim.com)