



BioStim[®] NMS²

Electrical Muscle Stimulator

Features

- ▶ Dual channel (4-lead) all digital Muscle Stim
- ▶ Three treatment modes (Cycled, Reciprocating, Constant)
- ▶ Choice of waveform (symmetrical or asymmetrical biphasic square wave)
- ▶ 8 Preset Programs
- ▶ 3 Customizable programs
- ▶ Fully adjustable rate, width, ramp, contraction and relaxation times
- ▶ Last setting recall & settings lock
- ▶ Large LCD display showing all parameters
- ▶ Compliance monitor
- ▶ Patient Lock
- ▶ Treatment timer
- ▶ Stores Custom Settings
- ▶ Removable plastic belt clip
- ▶ 5-year warranty

Technical Specifications

Channel: Dual

Amplitude: Continuously adjustable from 0 - 120 mA peak

Output Voltage: 0 - 60V peak to peak

Waveform: Symmetrical biphasic square, asymmetrical biphasic square

Pulse Rate: Adjustable 1 - 120 Hz

Pulse Width: Adjustable 50 - 400 μ S

Modes: Constant, Cycled, Reciprocating

Contraction Time: 0 - 99 seconds

Relaxation Time: 0 - 99 seconds

Ramp: Up/Down 0 - 10 seconds

Power Supply: 4 AA batteries, type LR6

Size: 9.9 x 6.5 x 3.3 cm

Weight: 226 grams

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 48" pin lead wires, 4 each AA alkaline batteries, carrying case, and instruction manual



Distributed by:

Medi-Stim, Inc., 217 Industrial Court
Wabasha, MN 55981 USA
Ph. 1-800-363-7846 Fax 1-651-565-2410
www.medi-stim.com sales@medi-stim.com