



CES≈ULTRA™

CRANIAL ELECTRICAL STIMULATION DEVICE

- ▶ Nonpharmacological alternative for the treatment and prevention of anxiety, depression and insomnia
- ▶ User-friendly: ease of application, single control setting
- ▶ Personal timer
- ▶ Efficient: automatic shut-off, low battery indicator
- ▶ Comfortable: choice of pre-gelled electrodes or ear clips
- ▶ Convenient: compact, portable, hand-held or placed in pocket
- ▶ Secure: rugged construction, flame-retardant, fitted carrying case

Technical Specifications:

Amplitude	Adjustable from 0 - 1.5 mA
Waveform	Square wave
Pulse Rate	100 Hz
Pulse Width	2 milliseconds (20% duty cycle)
Timer	30, 45 minutes, or continuous
Power Supply	9 volt alkaline battery
Size	107 x 65 x 25 mm
Weight	128 grams including battery

Distributed by

