

Available Exclusively from



Your Complete Source for Quality Electrotherapy Equipment and Supplies!

COMP STIM

FEATURES:

- Dual channel (4-lead) analog muscle stim
- 3 Treatment Modes (Alternating, Synchronous, or Constant)
- Adjustable frequency
- Adjustable contraction and relaxation times
- Adjustable ramp times
- 3 selectable pulse width settings (100, 200, or 300 μ S)
- Channel output indicator LED lights
- Locking front control panel
- Removable plastic belt clip
- 1-year warranty

Technical Specifications:

- ▶ Channel Dual, isolated between channels
- ▶ Amplitude Adjustable, 0 - 80 mA peak into 500 ohm load each channel
- ▶ Output Voltage Adjustable, 0 - 40 V peak into 500 ohm load each channel
- ▶ Waveform Asymmetrical biphasic square pulse
- ▶ Pulse Rate Adjustable, 2 - 150 Hz
- ▶ Pulse Width Adjustable, 30 - 260 μ S
- ▶ Modes Three (Normal, Burst, Modulation):
 - ▶ Burst (B) - Bursts occur twice per second. Pulse width adjustable, frequency fixed at 100 Hz
 - ▶ Modulation (M) - Pulse rate is automatically varied in a cyclic pattern over an interval of nominally 10 seconds (max. 150 Hz) Pulse rate decreases linearly over a period of 4 seconds from the control setting value to a value which is 40% less. The lower pulse rate will continue for 1 second, then increase linearly over a 4 second period to its original value. The original pulse rate will continue for 1 second. The cycle is then repeated.
- ▶ Timer 15 minutes, 30 minutes, or continuous
- ▶ Power Supply One 9-volt battery
- ▶ Size 9.5 x 6.5 x 2.35 cm
- ▶ Weight 115 grams with battery



Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43" pin lead wires, 1 each 9V heavy duty battery, carrying case, and instruction manual

Distributed by:



Medi-Stim, Inc.

217 Industrial Court, Wabasha, MN 55981 USA
Toll Free: 1-800-363-STIM | Fax 1-651-565-2410
www.medi-stim.com | sales@medi-stim.com