



# COMBO-STIM

## NEURO MUSCULAR ELECTRICAL STIMULATION UNIT

- + **2 Channels**
- + **2 Waveforms**
  - Symmetrical biphasic rectangular*
  - Alternating symmetrical biphasic rectangular*
- + **TENS FEATURES:**
  - 3 Programs (all fully adjustable): Constant, Burst, Modulated Rate & Width*
- + **EMS FEATURES:**
  - 5 Programs: 1 fixed preset and 4 adjustable presets*
  - 2 Modes: Synchronous, Alternating*
  - Adjustable Ramp UP and Ramp Down*
  - Channel 2 delay feature*
  - Stim On/Off Channel Indicators*
- + **Adjustable Pulse Rate & Pulse Width**
- + **Parameter and Program Locking**
- + **Fully digital with large LCD Display**
- + **2-year warranty**



### TECHNICAL SPECIFICATIONS

*Channel* Dual channel  
*Intensity* In 20 steps of 5mA, 5 to 100mA over 1000 ohm  
*Frequency* 1, 2, 3, 4, 5, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50, 60, 70, 80, 90, 110, 120, 130, 140, 150 Hz  
*Waveform* Symmetrical biphasic rectangular, Alternating symmetrical biphasic rectangular  
*Pulse Width* 50 - 400 $\mu$ s, in steps of 50 $\mu$ s  
*Ramp Time* 0 - 5s, in steps of 1s  
*On Time* 1 - 40s, in steps of 1s  
*Off Time* 0 - 40s, in steps of 1s  
*Treatment Timer* Continuous, 10, 20, 30, 45, 60, 90 minutes  
*Power Supply* 2 AA batteries

Distributed by:



217 Industrial Court  
Wabasha, MN 55981  
1-800-363-STIM | www.medi-stim.com  
Fax: 1-651-565-2410