

# ComforTrac

Distributed by:

**MEDI-STIM**.INC.

217 Industrial Court  
Wabasha, MN 55981  
1-800-363-STIM | www.medi-stim.com  
Fax: 1-651-565-2410

*Comfortable, simple & effective...*

## Lumbar Traction



The ComforTrac lumbar home traction unit has been carefully engineered to provide maximum comfort, efficacy and ease of use.




# ComforTrac

Distributed by:  
**MEDI-STIM** INC.  
217 Industrial Court  
Wabasha, MN 55981  
1-800-363-STIM | www.medi-stim.com  
Fax: 1-651-565-2410

*Our innovative design makes it possible for patients to go home with the same quality traction they receive during physical therapy.*

## Lumbar Traction



### SPECIFICATIONS:

**Maximum Force:** 200 lbs.

**Hand Pump:** Displays force in both pounds and kilograms

**Headrest Bed Material:** Memory foam

**Testing:** 5,000 cycles which equals using 3 times daily for 5 years

### SUPPORTING RESEARCH & REFERENCES:

Komori H, Shinomiya K, Nakai O, et. al., The Natural History of Herniated Nucleus Pulposus with Radiculopathy, Spine 21:225-229, 1996.

Onel D, Tuzlaci M. Sari H, et. al., Computed Tomographic Investigation of the Effect of Traction on Lumbar Disc Herniations, Spine 14:82-90, 1989

Matthews JA, Hickling J, Lumbar Traction: A Double-Blind Controlled Study for Sciatica. Rheumatology and Rehabilitation 14:222-225, 1975