

Electro-Mesh™ Conductive Garments

Electro-Mesh™ Conductive Electrode Garments contain the highest level of silver fibers compared to other brands. These top quality garments, made in USA, also contain patented Dacron™ fibers for optimum tensile strength, prolonged elasticity, and body heat retention (for increased circulation). As form-fitting garments, they provide stimulation throughout the entire fabric, offering greater stimulation coverage than conventional electrodes, and due to the silver, are ideal for patients who have skin reactions to pre-gelled electrodes. For best results, apply Electromist Conductive Spray to treatment area prior to wearing garment. *Note: Due to variations in garments, the stretch values listed are only guidelines. Patients at the high end of the stretch limit may want to try the next size up for comfort.*

Electro-Mesh™ Glove

- EMG2011** Extra Small (3.5" W x 7.0" L), size 5
- EMG2012** Small (4" W x 7.5" L), size 6
- EMG2013** Medium (4.5" W x 8.5" L), size 7
- EMG2014** Large (5" W x 9" L), size 8
- EMG2015** X-Large (5.5" W x 9.5" L), size 9



Electro-Mesh™ Sock

- EMS2112** Small (women's 6-9; men's 5-7)
- EMS2113** Medium (women's 10-11; men's 8-10)
- EMS2114** Large (women's 12-13; men's 11-12)
- EMS2115** X-Large (women's 14-15; men's 13-14)



Electro-Mesh™ Sleeve

- ESL2206** Small (6" L; 14" circumference over treatment area)
- ESL2208** Medium (8" L; 17" circumference over treatment area)
- ESL2210** Large (10" L; 20" circumference over treatment area)
- ESL2212** X-Large (12" L; 23" circumference over treatment area)



Electro-Mesh™ Wrap

- CEB2024** 2" W x 24" L Wrap
- CEB3036** 3" W x 36" L Wrap
- CEB4048** 4" W x 48" L Wrap



Electro-Mesh™ Silver Panels

Self-adhesive conductive garment fabric pads glued to a foam-backed panel. Peel off adhesive liner and apply pad directly to the inside of our conductive wraps or supports. Economical alternative to the standard Conductive Garment Electrodes. Apply conductive spray to surface before application for best results. Snap connector with pin adapter to accept snap or male 2mm pin from lead wire. Lasts approximately 2-3 months.

- EMSP002** Silver Panel Electrodes 2" x 2" (1 pair/pkg)
- EMSP003** Silver Panel Electrodes 3" x 4" (1 pair/pkg)
- EMSP008** Silver Panel Electrodes 8" x 4" (1 pair/pkg)



Care Vest Conductive Garment Electrode Vest

Allows patients to position electrodes, by themselves, in hard-to-reach areas such as the upper and mid-back regions. Constructed of soft, breathable materials for daily comfort. Hand or machine wash for easy care. Form-fitting, with full zipper front, and light-weight with pre-punched holes in the back panel for electrode connectors. Includes 1 pkg of 4 Stay-Put 2" x 2" velcro electrodes. Custom sizes available.

- CEV-S** Small chest: 34"-38" waist: 30-36"
- CEV-M** Medium chest: 38"-44" waist: 32-38"
- CEV-L** Large chest: 42"-48" waist: 36-40"
- CEV-XL** XL chest: 46"-52" waist: 38-42"
- CEV-2X** 2-XL chest: 50-60" waist: 40-46"



Stay-Put Electrodes

Stay-Put electrodes are pre-gelled self-adhesive reusable electrodes with a strong velcro hook backing to stay securely in place on the Conductive Garment Electrode Vest, or ThermoSkin® Wraps & Supports. Available in standard gel with white foam top or a soft blue gel with tan tricot knit top for patients with sensitive skin.

- SP5000V** 2" Round, Foam Pin Wire, Pregelled Electrodes w/ Velcro 4/pkg
- SP5050V** 2" x 2" Square, Foam Pin Wire, Pregelled Electrodes w/ Velcro 4/pkg
- SP5090V** 2" x 3.5" Rectangle, Foam Pin Wire, Pregelled Electrodes w/ Velcro 4/pkg
- SP7500V** 3" Round, Foam Pin Wire, Pregelled Electrodes w/ Velcro 4/pkg
- SS175V** 1.75" x 1.75" Square, Pin Wire, for Sensitive Skin Electrodes w/ Velcro 4/pkg
- SS375V** 1.75" x 3.75" Rectangle, Pin Wire, for Sensitive Skin Electrodes w/ Velcro 4/pkg



Conductive Garment Electrodes

Conductive Garment Electrode pads are made with highly conductive silver fabric fiber material. They offer very low impedance for comfortable stimulation. Great for patients with sensitivity to conventional gelled electrodes. Velcro hooks on the back side securely adhere to the inside surface of conductive wraps & supports. Foam padded to retain water or conductive spray. Snap connector with pin adapter to accept snap or male 2mm pin from lead wire. Lasts approximately 6-9 months. Hand wash and air dry.

- CGP83294** 2" x 3" Rectangle (1 pair/pkg)



ThermoSkin® Wraps and Supports

ThermoSkin® Conductive Wraps provide a means of positioning electrodes to specific body parts and keeping them in place, while offering compression and support. They are made of neoprene and Trioxon™, a patented material which not only captures and retains natural body heat (for increased circulation), but also, unlike other products, allows the skin to breathe and remain dry. Use velcro-backed Conductive Garment or Stay Put electrodes with these ThermoSkin® products.

Conductive Neck Wrap

- CSN82221** XS (11.5" - 12.75" neck)
- CSN83221** S (13" - 14.25" neck)
- CSN84221** M (14.5" - 15.75" neck)
- CSN85221** L (16" - 17.25" neck)
- CSN86221** XL (17.5" - 19" neck)
- CSN87221** XXL (19.25" - 21" neck)



Conductive Lumbar Support

- CBB82227** XS (23.5" - 27.25" waist)
- CBB83227** S (27.5" - 31.75" waist)
- CBB84227** M (32" - 35.5" waist)
- CBB85227** L (35.75" - 39.5" waist)
- CBB86227** XL (39.75" - 44" waist)
- CBB87227** XXL (44.25" - 48.5" waist)
- CBB88227** 3XL (48.75" - 53" waist)
- CBB89227** 4XL (53.25" - 57.75" waist)
- CBB80227** 5XL (58" - 62.25" waist)



Conductive Shoulder Wrap

- CSS82230** XS (33" - 35" chest)
- CSS83230** S (35.25" - 37.5" chest)
- CSS84230** M (37.75" - 40.25" chest)
- CSS85230** L (40.5" - 43" chest)
- CSS86230** XL (43.25" - 47.25" chest)
- CSS87230** XXL (47.5" - 51.5" chest)



Conductive Open Knee Wrap

- CKW82284** XS (10.25" - 11.75" knee)
- CKW83284** S (12" - 13.75" knee)
- CKW84284** M (14" - 15.75" knee)
- CKW85284** L (16" - 17.75" knee)
- CKW86284** XL (18" - 19.75" knee)
- CKW87284** XXL (20" - 21.75" knee)
- CKW88284** 3XL (22" - 22.75" knee)
- CKW89284** 4XL (24" - 25.75" knee)
- CKW80284** 5XL (26" - 27.75" knee)



Silver-Thera™ Conductive Garments

The Silver-Thera™ Garments are two Electro-Mesh™ electrodes sewn together into one piece for easier application, with a non-conductive material stitched in between to separate the two electrodes. A prosthetic sock electrode is also available for amputee patients suffering from phantom limb pain or to speed recovery from surgery. FDA-approved and Medicare reimbursable. Hand-washable. Sold individually.

Silver-Thera™ Glove

(device shown not included)

STG2412 Small (4" W x 7.5" L)
glove size 6

STG2413 Medium (4.5" W x 8.5" L)
glove size 7

STG2414 Large (5" W x 9" L)
glove size 8



Silver-Thera™ Stocking

(device shown not included)

STS2312 Small (women's size 6-9;
men's size 5-7)

STS2313 Medium (women's size 10-11;
men's size 8-10)

STS2314 Large (women's size 12-13;
men's size 11-12)

STS2315 X-Large (women's size 14-15;
men's size 13-14)

**Patients with larger calves may require being fitted with an Electro-Mesh™ sock and sleeve for better results.*



Silver-Thera™ Prosthetic Sock Electrode

(device shown not included)

STP5012 Regular Short (Fits Prosth. Sock 1-2) 10-14" L

STP5013 Regular Medium (Fits Prosth. Sock 1-2) 16-20" L

STP5014 Wide Medium (Fits Prosth. Sock 3+) 16-20" L



For best results, apply Electromist Conductive Spray to treatment area prior to wearing garment.


Back-Stim Wrap

The Back-Stim Wrap is a lumbar support belt with 4 pre-positioned conductive steel snaps on the back panel. Lead wires from the stimulator attach to the snaps on the outside of the belt. Current is conducted through the snaps to special double-sided pre-gelled electrode pads that are applied to the smooth surface of the snaps inside the belt. The pads conduct the current to the skin, and can be repositioned. Ideal for patients who have difficulty placing electrodes in the lumbar region by themselves, or difficulty keeping conventional electrodes in place. Constructed with a breathable back panel and made of washable neoprene and nylon. Available in Standard and Large sizes, with optional 12" belt extender. Includes 1 pkg of 4 Ultra-Stim #US2020 double sided gelled electrodes, and 4 pin-to-snap adapters.

USBB-S Standard (25"-42" waist)

USBB-L Large (40"-54" waist)

USBB-EXT 12" Wrap extender panel



Ultra-Stim® and Uni-Patch® Garment Electrodes


Conductive electrodes with adhesive gel on both sides, designed for use with Back-Stim and Ultra-Stim Wraps.

US2020 2" x 2" Square 4/pkg

US2040 2" x 4" Rectangle 2/pkg

PC90060 2" x 2" Square 4/pkg

PC91000 3.5" x 3.5" Square 4/pkg



Distributed by:



Conductive Garments



Electrotherapy – where you need it!

Conductive Garments are the newest method of applying electro-stimulation to targeted treatment sites of the body. Medi-Stim offers a full line of conductive garments that are available as socks, gloves, sleeves, wraps, and electrode pads, as well as form-fitted vests, lumbar belts, and orthopedic supports for other parts of the body. They can be used with all types of electro-therapy devices including TENS, muscle stim, high-volt galvanic, micro-current and interferential stimulators. Conductive Garments are covered by Medicare and most insurance companies. HCPC: E0731

**All of the products featured in this brochure are latex-free.*

Silver-Thera™ Garments

Back-Stim Wrap & Electrodes

Conductive Garments