

DURATENS

PERSONAL PAIN RELIEF SYSTEM

with high-tech design & protective holster

- Dual channel all-digital TENS
- Five treatment modes (Normal, Burst, Modulated Rate, Modulated Width 1, Modulated Width 2)
- Snap-on faux leather protective holster with metal belt clip
- Attractive high-tech design with metallic finish and LCD display
- Easy-to-use push-button controls for all settings
- Last setting recall & settings lock
- Compliance monitor
- Power ON/OFF button for easy shut-off
- LCD display shows all parameters on screen
- Treatment timer (30, 60 minutes or Continuous)
- 1 year warranty



Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43 p in lead wires, 1 each 9V heavy duty battery, carrying case, holster, and instruction manual



Technical Specifications

Channel Dual, isolated between channels
 Amplitude 0 - 80 mA peak (500 ohm load each channel)
 Output Voltage 0 - 40 V peak (500 ohm load each channel)
 Waveform Asymmetrical biphasic square pulse
 Pulse Rate Adjustable, 1 - 150 Hz
 Pulse Width Adjustable, 30 - 300 μ S
 Modes Five (Burst, Normal, Modulation, M1, M2)
Burst (B) Burst Rate adjustable 0.5 - 5 Hz, Pulse Width adjustable 30 - 300 μ S, Frequency fixed at 100 Hz
Normal (N) Pulse Rate adjustable 2 - 150 Hz, Pulse Width adjustable 30 - 300 μ S
Modulated Rate Pulse Rate decreases linearly over 5 seconds 40% of original setting, then increases over 5 seconds back to its original setting.
Modulated Pulse Width 1 decreases linearly over 5 seconds 40% of Width1: original setting, then increases over 5 seconds back to its original setting.
Modulation Pulse Width 2 decreases linearly over 5 seconds 70% of Width2: original setting, then increases over 5 seconds back to its original setting.
 Timer 30, 60 minutes, and continuous
 Power Supply One 9-volt battery
 Size 6.0 cm (L) x 9.3 cm (W) x 2.5 cm
 Weight 110 grams with battery

Distributed by:



217 Industrial Court
 Wabasha, MN 55981
 1-800-363-STIM | www.medi-stim.com
 Fax: 1-651-565-2410