

EMS-2C

DIRECT CURRENT POINT STIMULATOR

FEATURES

- Single channel (2-lead) analog direct current (DC) point stimulator
- Choice of Interrupted DC or Pulsed DC modes
- Hand-held trigger point or conventional electrode operation
- LED output and low battery indicators
- Cloth-covered treatment probe tip must be wet with water
- Ideal for facial and smaller peripheral muscles
- Reusable padded cloth dispersive grounding pad must be wet with water
- Push-button control to activate stimulation
- 2-Year Warranty

TECHNICAL SPECIFICATIONS

- Waveform: Monophasic
- Modes: Interrupted DC and Pulsed DC
 - Interrupted DC Amplitude 0-18 mA peak, Pulse Width 30 μS, Pulse Rate 1 Hz
 - Pulsed DC Amplitude 0-25 mA peak,
 Pulse Width 250 µS, Pulse Rate 80 Hz
- Power Supply: 9V battery
- Size: 140 x 108 x 41 mm
- Weight: 142 grams

Comes with 1 pair of reusable carbon rubber electrodes, 1 each grounding dispersive pad, 1 each 36 p in lead wire, 1 each 9-volt alkaline battery & instruction manual.

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA 1-800-363-STIM | Fax: 1-651-565-2410 | www.medi-stim.com



