

# ETPS 1000

## TENS/MENS Point Stimulator

### Features:

- Point finder and treatment probe in one tip
- Audible variable tone indicator for locating trigger points
- Remote extension probe for treating hard-to-reach areas
- Alternative for acupuncture therapy without needles
- Positive (tonify), Negative (disperse) or Bipolar (tonify & disperse) switch
- Trigger point Search and Stim button
- Adjustable skin sensitivity dial
- Battery test button
- LED output indicator light
- Adjustable wrist grounding strap
- 1-year warranty
- Comes with 1 each 9V alkaline battery, remote extension probe, wrist grounding strap, carrying case, and "Natural Health" Instruction & Protocol Booklet with illustrations and protocols for properly treating points for over 100 different pain disorders
- Leadwire reorder #TENSPR-WS (grounding wrist strap), #ETPS-RP (remote probe)

### Technical Specifications:

- Search Mode
- Output Parameters (skin resistance 500K ohms): Voltage 3.8 mv; Pulse Length 0.0909 microseconds; Frequency 5.518 KHz; Duty Cycle 51%
- Output Parameters (skin resistance 10K ohms): Voltage 71.6 mv; Pulse Length 0.0909 microseconds; Frequency 5.88 KHz; Duty Cycle 51%
- Stimulation Mode: Monophasic
- Output Parameters (resistance 500 ohms): Voltage .896 V; Pulse Length 213 milliseconds; Frequency 2.35 Hz; Duty Cycle 51%
- Output Parameters (resistance 10K ohms): Voltage 5.25 V; Pulse Length 213 milliseconds; Frequency 2.30 Hz; Duty Cycle 51%
- Stimulation Mode: Biphasic
- Output Parameters (resistance 500 ohms): Voltage 1.38 V; Pulse Length 213 milliseconds; Frequency 1.17 Hz; Duty Cycle 50%
- Output Parameters (resistance 10K ohms): Voltage 10.9 V; Pulse Length 213 milliseconds; Frequency 1.0 Hz; Duty Cycle 50%
- Maximum Pulse Width in Stimulation Mode 213 milliseconds for all polarity settings
- Maximum Intensity Setting 500 ohms resistance gives 85 microamps; 10K ohms resistance gives 600 microamps
- Minimum Intensity Setting 10K ohms resistance gives 300 microamps
- Power Source 9-volt battery

**MEDI-STIM** INC.

Medi-Stim, Inc.

217 Industrial Court, Wabasha, MN 55981 USA  
Ph. 1-800-363-STIM Fax 1-651-565-2410  
www.medi-stim.com sales@medi-stim.com

Distributed by:

