

# DYNATRON ibox™ Iontophoresis Unit

## Features:

- ▶ Dual channel digital iontophoresis drug delivery device
- ▶ Alternative to hypodermic needle injections of pain medication
- ▶ Independent channels to treat two areas with separate parameters
- ▶ Single phase (one drug) or dual phase (two drugs of opposing polarities delivered from the same iontophoresis electrode) from each channel
- ▶ Positive(+), negative(-), and dual polarity (+/-) options
- ▶ Bright OLED display shows channels, intensity, dosage and treatment time
- ▶ Screen shows error messages for troubleshooting assistance
- ▶ Energy efficient: delivers max. 3,360 mA-min. with 2 AA batteries
- ▶ Custom presets & defaults for each channel
- ▶ Detects electrode / lead errors
- ▶ 2 year warranty
- ▶ Comes with 2 pairs of 36" snap squeeze-clip color-coded lead wires, 2 AA alkaline batteries, carrying case, and instruction manual

\*Electrodes sold separately.

## Technical Specifications:

Output Jacks 2  
 Current (Intensity) 0.5 – 4.0 mA  
 Current Accuracy +/- 10%  
 Current Precision 0.1 mA  
 Voltage at Electrodes +/- 0-46 volts  
 Electrode Cable Jacks Custom 2-pin polarized  
 Maximum Total Deliverable Dose 3360 mA-min.  
 (using high quality Alkaline batteries)  
 Single Treatment Deliverable Dose 1-160 mA-min.  
 Number of 40 mA-min. Treatments 89  
 Treatment Time 15 sec. - 100 minutes  
 Timing Accuracy +/- 1 second  
 Certificates IEC 60601-1 (CE) and CSA/NRTL  
 Display Organic light emitting diodes (OLED)  
 Size 4.5" x 2.7" x 1.25"  
 Weight 3.6 oz (w/o battery or electrode cables)  
 Power Supply 2AA batteries



Use with Medi-Stim IO Electrodes!



# IONTO15  
Small (1.5cc/4.0mA)  
2.5" Round

# IONTO25  
Medium (2.5cc/4.0mA)  
3" Round

# IONTO40  
Large (4.0cc/4.0mA)  
3.25" Round

# IONTO20B  
Butterfly (2.0cc/4.0mA)  
3"x 4" Rectangle

## Medi-Stim IO Electrode Features:

- ▶ Buffered return electrode with excellent current dispersion and lower impedance for effective delivery of solution
- ▶ Larger surface area of return electrode provides greater comfort during treatment
- ▶ Unique buffered gel maintains optimum pH levels
- ▶ Soft, skin friendly, and conformable materials for buffered active electrode; offers excellent adhesion to minimize leakage
- ▶ Round shape design of electrodes provides improved conformity and adhesion to the treatment site
- ▶ User-friendly positioning grip tabs for easy application and removal
- ▶ Return and active electrodes packed together in long-life Mylar foil package
- ▶ Standard snap connector
- ▶ 12 kits/box; each kit contains one active electrode and one buffered electrode

Distributed by:



217 Industrial Court  
 Wabasha, MN 55981  
 1-800-363-STIM | www.medi-stim.com  
 Fax: 1-651-565-2410

**Caution:** Federal law restricts this device to sale by or on the order of a physician.