



# Impulse HVG

High Volt Pulsed Galvanic Stimulator



## Features:

- ▶ Single channel (2-lead) digital high volt pulsed galvanic stimulator (HVPG)
- ▶ Manual or one of 3 preset protocols (Chronic pain, Edema reduction, or 48-hour Post-op)
- ▶ Tri-lead cable with 1 dispersive & 2 active leads
- ▶ Alternating or continuous treatment mode switch
- ▶ Negative or positive polarity switch
- ▶ Adjustable frequency (1-120 Hz)
- ▶ Energy-efficient micro-processor technology
- ▶ Compliance monitor
- ▶ Removable plastic belt clip
- ▶ Limited lifetime warranty

## Technical Specifications:

**Channel:** Single, 1 dispersive and 2 active leads

**Output Voltage:** Min. 200V - 300V

**Waveform:** Monophasic pair of spiked waves

**Pulse Rate:** 1-120 pulses per second

**Power Supply:** One 9V alkaline battery

**Size:** 11.4 x 6.3 x 1.9 cm

**Weight:** 90.7 grams (without battery)



Comes with 1 each 4x7 s self-adhesive dispersive electrode, 4 self-adhering 1.5 x 2.5 ov al active electrodes, 36 T ri-lead cable with 1 dispersive lead and 2 active leads, 1 each 9V alkaline battery, carrying case, and instruction manual.



Your Complete Source for Quality Electrotherapy Equipment and Supplies!

217 Industrial Court, Wabasha, MN 55981 USA  
Ph. 1-800-363-STIM Fax 1-651-565-2410  
www.medi-stim.com sales@medi-stim.com