


MEDI-STIM INC.

MH8000

Combination TENS/NMES Unit

Features:

- ▶ Dual channel (4-lead) digital combo TENS and muscle stim (EMS)
- ▶ 5 TENS modes (Burst, Normal, Modulated Width, Modulated Rate, and Modulated Rate & Width)
- ▶ 2 EMS modes (Simultaneous, Alternating)
- ▶ All digital push-button controls
- ▶ Treatment timer (5-90 minutes or Continuous)
- ▶ Last setting recall & settings lock
- ▶ Stimulation output indicators
- ▶ Compliance monitor
- ▶ Low battery power indicator
- ▶ Hi-tech design with metallic finish w/ large LCD display
- ▶ Removable plastic belt clip
- ▶ 1-year warranty
- ▶ Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43" pin lead wires, 1 each 9V heavy duty battery, carrying case, and manual
- ▶ Leadwire reorder #LW46SPLE



Technical Specifications:

Channel: Dual, isolated between channels.

Amplitude: Adjustable 0 - 80 mA, Max output 100 mA peak to peak (21.2 mA rms) into 500 ohm load each channel.

Output Voltage: Adjustable 0 - 50 V, Max output 50 V peak to peak (10.6V rms) into 500 ohm load each channel.

Waveform: Asymmetrical rectangular biphasic pulse.

Pulse Rate: Adjustable, 2 - 20 Hz in 1 Hz/step and 20-150 Hz in 5 Hz/step.

Pulse Width: Adjustable, 50µS - 300µS, 10µS/step.

5 TENS Modes: Burst, Normal, Modulated Width, Modulated Rate, and Modulated Rate & Width.

2 NMES Modes: Simultaneous or Alternating.

NMES Contraction Time: Adjustable, 2 - 90 seconds, 1 Sec./ step.

NMES Relaxation Time: Adjustable, 2 - 90 seconds, 1 Sec./ step.

NMES Ramp Time: Adjustable, 1 - 8 seconds, 1 Sec./ step.

Timer: Adjustable, 5 - 90 minutes or Continuous.

Power Supply: One 9-volt battery.

Size: 13.6(L)x7(W)x2.65(H)(cm).

Weight: 162 grams (battery included)



Distributed by:

MEDI-STIM INC.

217 Industrial Court
Wabasha, MN 55981

1-800-363-STIM | www.medi-stim.com
Fax: 1-651-565-2410

Medi-Stim, Inc., 217 Industrial Court, Wabasha, MN 55981 USA
Ph. 1-800-363-STIM Fax 1-651-565-2410
www.medi-stim.com sales@medi-stim.com

Caution: Federal law restricts this device to sale by or on the order of a physician.