

# MS-TENS+

With AC Adaptor!

- Dual channel (4-lead) analog TENS
- Three treatment modes (Normal, Burst, Modulated Width)
- Treatment timer (30, 60 minutes or Continuous)
- Easy-to-use dials and controls
- Protective slide-down control panel
- Recessed amplitude dials
- Powered by 9V battery or AC adaptor
- Removable plastic belt clip
- 1-year warranty
- Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43" pin lead wires, 1 each 9V heavy duty battery, zippered case pouch, and instruction manual
- Leadwire reorder #LW46SPL



## Technical Specifications:

**Channel** Dual, isolated between channels

**Amplitude** 0 - 80 mA

**Output Voltage** 0-45V peak to peak into 500 ohm load each channel

**Waveform** Asymmetrical Biphasic Square Wave

**Pulse Rate** 2 - 150 Hz

**Pulse Width** 40 - 260  $\mu$ S

**Modes** Burst, Normal and Modulated Width

**Treatment Timer** 30 or 60 minutes or Continuous

**Power Source** 9 volt battery or AC adaptor

**Size** 9.4 x 6.4 x 3 cm

**Weight** 136 grams with battery



Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43" pin lead wires, AC Adaptor, 1 each 9V heavy duty battery, carrying case, & instruction manual

Distributed by:



217 Industrial Court  
Wabasha, MN 55981

1-800-363-STIM | [www.medi-stim.com](http://www.medi-stim.com)

Fax: 1-651-565-2410