



# MICRO 400

## MENS/TENS Unit

### Features:

- *Dual channel (4-lead) analog microcurrent*
- *Selectable frequencies (.3, 9.1, 30, or 292 Hz)*
- *Choice of micro-ampere or milli-ampere (TENS) currents*
- *Square or sloped waveform*
- *Adjustable amplitude*
- *Bipolar [(+) and (-)] direct current waveform*
- *LED output indicator light*
- *Durable metal belt clip*
- *1-year warranty*



### Technical Specifications:

<b>Channel:</b>	Dual, isolated between channels
<b>Amplitude:</b>	Adjustable, 0 - 8 mA
<b>Waveform:</b>	Square or slope
<b>Carrier Frequency:</b>	15000 Hz
<b>Modulating Frequency:</b>	0.3Hz, 9.1Hz, 30Hz, 292Hz
<b>Power Supply:</b>	9-volt alkaline battery
<b>Size:</b>	2.7 x 6.2 x 9.5 cm
<b>Weight:</b>	150 grams including battery

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43 pin lead wires, 1 each 9V heavy duty battery, carrying case, and instruction manual.

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA  
1-800-363-STIM | Fax: 1-651-565-2410 | [www.medi-stim.com](http://www.medi-stim.com)