


# mini TENS/EMS



## features

- Single channel (2-lead) all digital combo TENS and EMS
- Snap-on "butterfly" electrode for wireless use
- 36" pin cable can be snapped on main unit to use with conventional electrodes
- 5 Preset TENS programs with timed treatments
- 5 Preset EMS programs with timed treatments
- Compact and lightweight design with watch battery
- Allows patient to discreetly wear under clothing
- Easy to use and set up
- 1-year warranty

## technical specifications

- Channel Single channel output
- Density Change of 16 steps pulses amplitude
- Output Voltage 0 ~ 90 mA peak to peak (load 500 ohm)
- Waveform Dual phase square wave
- Display LCD display with timer/mode/density
- Modes: Part A - P1~P5 are 5 TENS program  
Part B - P6~P0 are 5 EMS program; Rate: 2 Hz~86 Hz;  
Contraction (on time) 2 sec. or 4 sec.; Relaxation  
(off time) 2 sec. or 4 sec.
- Timer 30 minutes backward to auto turn off
- Power Supply 3V CR2025 battery
- Size 5.3 x 3 x 1 cm
- Weight Less than 30 grams including battery

Comes with 1 each reusable self-adhesive 2" x 4" butterfly mini-snap electrodes, 1 each 36" mini-snap lead wire, 1 each 3V lithium battery, plastic storage box, and instruction manual.

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA  
1-800-363-7846 | Fax: 1-651-565-2410 | www.medi-stim.com