

Myo-Care EMS



- Dual channel (4 lead) digital muscle stim with LCD display
- Synchronous mode (both channels contract and relax at the same time)
- Treatment timer prevents over-stimulation
- Function indicator
- Adjustable rate & width, ramp, & contraction/relaxation
- Hinged protective cap for intensity dials
- Easy-to-use push button controls
- Last setting recall
- Compliance monitor
- Metallic finish w/ durable wire metal clip
- Three year warranty
- Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43" pin lead wires, one 9V heavy duty battery, carrying case, and instruction manual



Technical Specifications

Channel Dual, isolated between channels
Amplitude 0 - 90 mA peak to peak (500 ohm load)
Output Voltage 0 - 40 V peak (500 ohm load)
Waveform Asymmetrical biphasic square pulse
Pulse Rate Adjustable, 1Hz ~ 150 Hz, 1 Hz/step
Pulse Width Adjustable, 30 μ S ~ 260 μ S, 10 μ S/step
Timer Adjustable, 1~ 60 minutes, 1min./step
Contraction Time Adjustable, 1~ 30 sec., 1 sec./step
Relaxation Time Adjustable, 1~ 30 sec., 1 sec./step
Ramp Up/Down Adjustable, 1~ 6 sec., 1 sec./step
Power Supply One 9-volt battery
Size 10 x 6 x 2.1 cm
Weight 120 grams including battery

Distributed by:



217 Industrial Court
Wabasha, MN 55981
1-800-363-STIM | www.medi-stim.com
Fax: 1-651-565-2410