



*No Prescription Required  
FDA Cleared!*

# PainAElly

**Non-Prescription TENS/NMS Device**

## FEATURES:

- ▶ Dual channel digital combination TENS & EMS device
- ▶ FDA approved for Over-the-Counter use; no Rx required
- ▶ 7 Preset TENS programs
- ▶ 6 Preset EMS programs
- ▶ Preset program treatment times (TENS: 1 hour ; EMS: 15 minutes)
- ▶ Open electrode detection function
- ▶ Low battery indicator
- ▶ 1 year warranty

## TECHNICAL SPECS:

Channels: Dual, isolated

Waveform: Asymmetrical, rectangular bi-phasic with zero DC current.

Current: Constant current

Amplitude: 0-80 mA

Dimensions: 108 x 62 x 23 mm.

Weight: 0.07KG without battery, 0.1KG with battery.



## TENS

- Program selectable pulse width: 100  $\mu$ S - 200  $\mu$ S
- Program pulse rate selection: In the continuous mode 10  $\square$  150 Hz
- Mode: Continuous, Burst, Modulation.
- Burst mode: Bursts of 9 pulses 150 Hz at 200  $\mu$ S, repeating twice every second.
- Modulation mode: 6-second cycle of concurrent pulse width modulation and rate modulation.
- Time duration of programs: 1 hour.

## NMS

- Program selectable pulse width: 200  $\mu$ S - 250  $\mu$ S
- Program pulse rate selection: in the continuous mode 12  $\square$  35 Hz
- Time duration of programs: 15 minutes

System includes: Stimulator, 1 set of 40" leadwires, 1 garment belt, 1 pkg of 4 reusable snap electrodes, 1 9V alkaline battery, instruction manual and carrying case.

