



SPORTX

COMBINATION PULSED GALVANIC
AND TENS STIMULATOR

▶ **Two Isolated Channels:**

Upper Channel: TENS (AC) to treat pain

Lower Channel: Pulsed Galvanic (DC) to reduce edema or increase circulation

▶ **Continuous or Preset program for ease of use**

▶ **Unique Triphasic Waveform for more comfort (AC Square Wave is added to the DC pulse)**

▶ **Push-Button Controls for easy use**

▶ **Compact, Portable and Lightweight**

▶ **Used by Professional Athletes**

▶ **3-Year Warranty**

Technical Specifications:

Channel: Dual, isolated between channels

Amplitude: Adjustable, 0 - 100 mA

Waveform: Biphasic, symmetrical square, interphase interval is 1800 μ S. Unbalanced, symmetrical triphasic square. Positive phase to negative phase interval is 1800 μ S, negative phase to negative phase is 900 μ S. Waveform period is 4500 μ S.

Pulse Rate: Biphasic waveform is 278 Hz, Triphasic waveform is 222 Hz.

Pulse Width: 60 μ S, fixed

Mode On second push of timer button, intensity ramps up to 90% of displayed level. After 31.5 minutes, stimulation ramps down 50% and returns to standby mode.

Programmed Time: Total treatment time of 39.5 minutes

Standby Mode: On first push of timer button, unit shuts off output and enters a standby mode. The unit will shut off after 8 hours of no use.

Maximum Charge per Pulse: 6 Microcoulombs

Power Supply: One 9-volt rechargeable or alkaline battery

Size: 27.1 x 13.5 x 5.9 cm

Weight: 124 grams including battery

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 36" pin lead wires, 2 each 9V alkaline batteries, zippered carrying case, and instruction manual



Distributed by:



217 Industrial Court Wabasha, MN 55981 USA
1-800-363-7846 | Fax: 1-651-565-2410 | www.medi-stim.com