



TENS PRO 900

Features:

- Hand-held acupuncture/trigger point stimulator
- Point finder and treatment probe in one tip
- Audible variable tone indicator for locating trigger points
- Remote extension probe for treating hard-to-reach areas
- Alternative for acupuncture therapy without needles
- Positive (tonify) or Negative (disperse) switch
- Trigger point Search and Stim button
- Adjustable skin sensitivity dial
- Battery test button
- LED output indicator light
- Adjustable wrist grounding strap
- 1-year warranty



Technical Specifications:

Waveform: Square waveform pulse
User-Selected Parameters
Tonify/Disperse Switch
Sensitivity Control: For optimizing detection of stimulation points by varying applied voltage.
Search/Stim Button
Intensity Control (in battery compartment): Reduces output in stimulation mode for the comfort of the patient (rarely adjusted).
Polarity: Positive or negative.
Search Mode Output Parameters (at a skin impedance of 100K ohms): Voltage 0.15 - 4V; Pulse Width 0.3-3.6 msec; Frequency 200 Hz - 2 KHz
Stim Mode Output Parameters (at a skin impedance of 100K ohms): Voltage 0-8.8V; Frequency 2.5 Hz
Power Supply: One 9-volt battery
Size: 18.8 x 6 x 2.4 cm
Weight: 122 grams

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA
1-800-363-STIM | Fax: 1-651-565-2410 | www.medi-stim.com

Comes with 1 each 9V alkaline battery, remote extension probe, wrist grounding strap, carrying case, and a very comprehensive Natural Health Instruction & Protocol Booklet with illustrations and protocols for properly treating points for over 100 different pain disorders.