



TS-1211

- ▶ Dual channel (4-lead) digital TENS
- ▶ Five treatment modes (Normal, Burst, Modulated Rate & Width, Strength Duration (SD) and Bi-Pulse)
- ▶ Treatment Timer
- ▶ Hinged protective cap for intensity dials
- ▶ Easy-to-use push-button controls
- ▶ Low-battery indicator
- ▶ Removable plastic belt clip
- ▶ 1-year warranty

TECHNICAL SPECIFICATIONS

Channel: Dual, isolated between channels

Amplitude: Adjustable 0 - 80 mA peak into 500 ohm load each channel

Output Voltage: 0 - 40V (load 500 ohm)

Waveform: Asymmetrical biphasic square pulse

Pulse Rate: Adjustable 2-150 pulses per sec. (Hz)

Pulse Width: Adjustable 50-300 μ S

Modes: Normal, Burst, Modulated Rate & Width, Strength Duration, Bi-Pulse

Patient Compliance Timer: Records the last 60 individual treatment sessions, with up to 999 minutes per session, and can record the whole treatment time for all the cumulative sessions.

Treatment Timer: Adjustable 5 - 90 minutes in 5 minute increments, or continuous

Power Supply: One 9-volt battery (heavy duty, alkaline, or rechargeable)

Size: 10 x 5.9 x 2.4 cm

Weight: 130 grams with battery



Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43" pin lead wires, 1 each 9V heavy duty battery, carrying case and instruction manual.

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA
1-800-363-7846 | Fax: 1-651-565-2410 | www.medi-stim.com