

ULTIMA 20

ULTIMA 20 TENS UNIT FEATURES:

- ▶ Dual channel (4-lead) digital TENS
- ▶ Three selectable waveforms (symmetrical biphasic, asymmetric biphasic, or monophasic)
- ▶ 20 treatment modes (8 manual modes and 12 preset modes)
- ▶ Last setting recall & settings lock
- ▶ Large LCD display screen with blue backlight
Treatment timer (15, 30, 45, 60, 90 minutes, or Continuous)
- ▶ Plastic belt clip
- ▶ 1-year warranty



TECHNICAL SPECIFICATIONS:

Channel: Dual

Amplitude Maximum: 130 mA peak value.

Output Voltage: 0 - 73 V across 500 ohm load.

Pulse Rate: 1Hz - 150Hz

Waveform: Symmetrical Biphasic rectangular, Monophasic rectangular, Asymmetrical biphasic rectangular.

Pulse Width: Adjustable, 50 - 250 μ S

Modes: 8 manual modes, 12 preset modes.

Timer: 15, 30, 45, 60, or 90 minutes, or continuous

Power Supply: 2 AA alkaline batteries

Size: 10.5 x 5.8 x 2.8cm

Weight: 136 grams with batteries

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43" pin lead wires, 2 AA alkaline batteries, carrying case and instruction manual. Leadwire reorder #LW46SPLE.



Distributed by:

Medi-Stim, Inc.

217 Industrial Court, Wabasha, MN 55981 USA

Ph. 1-800-363-STIM Fax 1-651-565-2410

www.medi-stim.com sales@medi-stim.com