



Ultima 11

DIGITAL TENS STIMULATOR

Features

- ▶ Dual channel (4-lead) all digital TENS
- ▶ Three selectable waveforms (asymmetrical or symmetrical biphasic, or monophasic rectangular)
- ▶ 11 Treatment modes (Constant 1, 2, or 3, Burst 1, 2, or 3, Modulated Width, Modulated Rate, Modulated Rate & Width, Chronic Preset, or Acute-Preset)
- ▶ Large LCD display with larger text
- ▶ Last Setting Recall & Settings Lock
- ▶ Treatment Timer (15, 30, 45, 60, 90 minutes, or Continuous)
- ▶ Compliance Monitor
- ▶ Removable plastic belt clip
- ▶ 1-year warranty

Technical Specifications

Channel: Dual

Amplitude: Maximum 130 mA peak across 500 ohm load

Output Voltage: 0 - 73 V across 500 ohm load

Waveform: Selectable, Asymmetrical or Symmetrical Biphasic, or Monophasic Rectangular

Pulse Rate: Adjustable, 1 - 150 Hz

Pulse Width: Adjustable, 50 - 250 μ S

Modes: Constant 1, 2, or 3, Burst 1, 2, or 3, Modulated Width, Modulated Rate, Modulated Rate & Width, Chronic (preset) and Acute (preset)

Timer: 15, 30, 45, 60, or 90 minutes, or continuous

Power Supply: 2 AA alkaline batteries

Size: 13.8 x 6.8 x 2.8 cm

Weight: 150.3 grams with batteries



Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA
1-800-363-STIM | Fax: 1-651-565-2410 | www.medi-stim.com

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43 p in lead wires, 2 each AA alkaline batteries, compact carrying case, and instruction manual