



# Ultima 11

DIGITAL TENS STIMULATOR

## Features

- ▶ *Dual channel (4-lead) all digital TENS*
- ▶ *Three selectable waveforms (asymmetrical or symmetrical biphasic, or monophasic rectangular)*
- ▶ *11 Treatment modes (Constant 1, 2, or 3, Burst 1, 2, or 3, Modulated Width, Modulated Rate, Modulated Rate & Width, Chronic Preset, or Acute-Preset)*
- ▶ *Large LCD display with larger text*
- ▶ *Last Setting Recall & Settings Lock*
- ▶ *Treatment Timer (15, 30, 45, 60, 90 minutes, or Continuous)*
- ▶ *Compliance Monitor*
- ▶ *Removable plastic belt clip*
- ▶ *1-year warranty*

## Technical Specifications

**Channel:** Dual

**Amplitude:** Maximum 130 mA peak across 500 ohm load

**Output Voltage:** 0 - 73 V across 500 ohm load

**Waveform:** Selectable, Asymmetrical or Symmetrical Biphasic, or Monophasic Rectangular

**Pulse Rate:** Adjustable, 1 - 150 Hz

**Pulse Width:** Adjustable, 50 - 250  $\mu$ S


**Modes:** Constant 1, 2, or 3, Burst 1, 2, or 3, Modulated Width, Modulated Rate, Modulated Rate & Width, Chronic (preset) and Acute (preset)

**Timer:** 15, 30, 45, 60, or 90 minutes, or continuous

**Power Supply:** 2 AA alkaline batteries

**Size:** 13.8 x 6.8 x 2.8 cm

**Weight:** 150.3 grams with batteries



Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA  
1-800-363-STIM | Fax: 1-651-565-2410 | www.medi-stim.com

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43 p in lead wires, 2 each AA alkaline batteries, compact carrying case, and instruction manual