



Ultima 5

DIGITAL TENS STIMULATOR

Features

- ▶ *Dual Channel (4-lead) Digital TENS*
- ▶ *Two Selectable Waveforms: (Symmetrical biphasic or Monophasic)*
- ▶ *Five Treatment Modes (Constant, Burst I, Burst II, Modulated Width, Modulated Rate & Width)*
- ▶ *Large LCD display showing all parameters*
- ▶ *Last Setting Recall & Settings Lock*
- ▶ *Treatment Timer (15, 30, 45, 60, 90 minutes, or Continuous)*
- ▶ *Removable plastic belt clip*
- ▶ *1-year warranty*

Technical Specifications

Channel: Dual
Amplitude Maximum: 150mA (peak value) across 500 ohm load
Output Voltage: 0-73 V across 500 ohm load
Waveform: Symmetrical Biphasic, or Monophasic
Pulse Rate: Adjustable, 1 - 150 Hz
Pulse Width: Adjustable, 50 - 250 μ S in steps of 10 μ S Treatment
Timer: 15, 30, 45, 60, or 90 minutes, or Continuous
Treatment Modes: Constant, Burst I, Burst II, Modulated Width, Modulated Rate & Width
Power Supply: 2 AA alkaline batteries 1.5 V
Dimensions: 10 x 6.5 x 2.8 cm
Weight: 127.5 grams (including batteries)



Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA
1-800-363-STIM | Fax: 1-651-565-2410 | www.medi-stim.com



Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43 pin lead wires, 1 each 9V heavy duty battery, carrying case, and instruction manual