

NO PRESCRIPTION REQUIRED!



Pain Warrior®

TENS/EMS Combo Stimulator



Features

- Dual Channel Digital Combination TENS/EMS Device
- FDA cleared; No Prescription Required
- 7 Preset TENS Programs
- 6 Preset EMS Programs
- Open Electrode Detection Function
- Low Battery Indicator
- 2 Year Warranty
- Leadwire Reorder: # LWSTP36

	Program	Description	Frequency	Width	On Time	Off Time	Time
TENS	1	CON	80	200			1 hr
	2	CON	50	100			1 hr
	3	CON	10	175			1 hr
	4	BST	150	200			1 hr
	5	BST	150	175			1 hr
	6	MOD	100/65	200/100			1 hr
	7	MOD	65/100	200/100			1 hr
EMS	8	NMS	12	200	5	5	15 min
	9	NMS	35	250	8	8	15 min
	10	NMS	12	200	5	10	15 min
	11	NMS	35	200	6	12	15 min
	12	NMS	12	250	5	15	15 min
	13	NMS	35	200	6	18	15 min



Technical Specifications

Channels: Dual, Isolated

Amplitude: 0-80 mA

Waveform: Asymmetrical, Rectangular bi-phasic, with zero DC current

Power Supply: One 9 Volt Battery

Size: 10.5 (L) x 6.5 (W) x 2.25 (D) cm

Weight: 116 grams (with battery)

System Includes

- 1 Pain Warrior Stimulator
- 1 set of 40" leadwires
- 1 garment belt
- 1 pkg of 4 reusable electrodes
- 1 9 Volt battery
- Protocol Index Cards
- Instruction Manual
- Carrying Case

Distributed by: