



# E-WAVE

## FEATURES

- ▶ Dual channel (4-lead) all digital muscle stim/FES
- ▶ Two treatment modes (simultaneous or alternating)
- ▶ Choice of waveform (symmetrical biphasic or pulsed direct current)
- ▶ Pulsed direct current waveform for edema reduction or increase circulation
- ▶ Fully adjustable rate, width, ramp, contraction and relaxation times
- ▶ External remote trigger switch input
- ▶ Last setting recall & settings lock
- ▶ Compliance Monitor
- ▶ Treatment Timer
- ▶ Removable plastic belt clip
- ▶ 3-Year warranty

## TECHNICAL SPECIFICATIONS

- Channel:** Dual channel  
**Amplitude:** Adjustable, 0 - 100 mA  
**Waveform:** Symmetrical Biphasic AC or Pulsed DC  
**Pulse Rate:** Adjustable, 4 - 100 Hz  
**Pulse Width:** Adjustable, 50 - 300  $\mu$ s  
**Modes:** Select either Simultaneous stimulation of both channels or alternating, stimulation 1 channel at a time  
**Timer:** Continuous operation or selectable 10-100 minutes in steps of 10 minutes  
**Ramp:** Up/Down 0.1 - 6 seconds  
**Power Supply:** 9-volt battery  
**Size:** 6.35 x 14 x 2.5 cm  
**Weight:** 227 grams including battery

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA  
1-800-363-7846 | Fax: 1-651-565-2410 | www.medi-stim.com



Safe & easy to use  
State-of-the-art  
technology!  
Powerful & comfortable  
stimulation.



Comes with 2 pkg of 4 reusable self-adhering electrodes, 1 set 40" pin lead wires, 1 each 9V alkaline battery, zippered carrying case and instruction manual. \*Available hand switch sold separately.