



E-WAVE

FEATURES

- ▶ Dual channel (4-lead) all digital muscle stim/FES
- ▶ Two treatment modes (simultaneous or alternating)
- ▶ Choice of waveform (symmetrical biphasic or pulsed direct current)
- ▶ Pulsed direct current waveform for edema reduction or increase circulation
- ▶ Fully adjustable rate, width, ramp, contraction and relaxation times
- ▶ External remote trigger switch input
- ▶ Last setting recall & settings lock
- ▶ Compliance Monitor
- ▶ Treatment Timer
- ▶ Removable plastic belt clip
- ▶ 3-Year warranty

TECHNICAL SPECIFICATIONS

- Channel:** Dual channel
Amplitude: Adjustable, 0 - 100 mA
Waveform: Symmetrical Biphasic AC or Pulsed DC
Pulse Rate: Adjustable, 4 - 100 Hz
Pulse Width: Adjustable, 50 - 300 μ s
Modes: Select either Simultaneous stimulation of both channels or alternating, stimulation 1 channel at a time
Timer: Continuous operation or selectable 10-100 minutes in steps of 10 minutes
Ramp: Up/Down 0.1 - 6 seconds
Power Supply: 9-volt battery
Size: 6.35 x 14 x 2.5 cm
Weight: 227 grams including battery

Distributed by:



Medi-Stim, Inc., 217 Industrial Court Wabasha, MN 55981 USA
1-800-363-7846 | Fax: 1-651-565-2410 | www.medi-stim.com



Safe & easy to use
State-of-the-art
technology!
Powerful & comfortable
stimulation.



Comes with 2 pkg of 4 reusable self-adhering electrodes, 1 set 40" pin lead wires, 1 each 9V alkaline battery, zippered carrying case and instruction manual. *Available hand switch sold separately.