



Distributed by:
MEDI-STIM INC.
217 Industrial Court
Wabasha, MN 55981
1-800-363-STIM | www.medi-stim.com
Fax: 1-651-565-2410

Comfortable, simple & effective...

Cervical Traction




The ComforTrac cervical home traction unit has been carefully engineered to provide maximum comfort, efficacy and ease of use.

Our innovative design makes it possible for patients to go home with the same quality traction they receive during physical therapy.


Cervical Traction

Fully Adjustable Angle of Incline
without the need to order extra parts



Hand Pump
our patented hand pump makes it easy to increase or decrease traction force

Quick Release Button
push to release the traction force



Memory Pillow
conforms comfortably to the head

Independently Adjustable Neck Wedges
provide a custom fit

SPECIFICATIONS:

Maximum Force: 50 lbs.

Angles of Incline: 10, 15 and 20 degrees (with no additional parts or attachments)

Hand Pump: Displays force in both pounds and kilograms

Neck Wedges: Independently adjustable

Headrest Material: Memory foam

Testing: 15,000 cycles which equals using 3 times daily for 5 years

SUPPORTING RESEARCH & REFERENCES:

Spine, 1996 Aug 15;21(16):1877-83 "Nonoperative Management of Herniated Cervical Intervertebral Disc with Radiculopathy." PMID: 8875719

Spine, 1985 Dec 10(10):867-71. "Cervical Spine Disorders. A Comparison of Three Types of Traction." PMID:3914085

J Manipulative Physiol Ther. 2002 Mar-Apr;25(3):188-92. "Intermittent Cervical Traction for Cervical Radiculopathy Caused by Large-Volume Herniated Discs." PMID: 11986581

Bangladesh Med Res Counc Bull. 2002 Aug;28(2):61-9. "Effects of Cervical Traction and Exercise Therapy in Cervical Spondylosis." PMID 12825763