


**THE NEW, SIMPLE & USER FRIENDLY**

# EZ-TENZ

**CHOOSE FROM 9 PRESET PROGRAMS**

**FOR CHRONIC/ACUTE/BURST PAIN**

**OR 6 BODY-SPECIFIC PRESETS**



## PROGRAMS

1. Chronic Pain
2. Acute Pain
3. Burst Mode
4. Neck
5. Shoulder
6. Lower Back/Hip
7. Hand/Elbow
8. Leg/Knee
9. Ankle/Foot

## FEATURES

- Dual channel (4-lead) digital TENS
- 9 Preset Programs defined for Body Regions
- Back Lit LCD Display
- Selectable Treatment Timer (15, 30, 60 min. or Continuous)
- Last Setting Recall
- Low Battery power indicator
- Hinged protective cap for intensity dials
- Removable Belt Clip
- 3 year warranty

### SYSTEM INCLUDES:

1 Stimulator, 4 Reusable Electrodes,  
1 Set 43" Lead Wires, 1 9V Battery,  
Carrying Case & Instruction Manual



217 Industrial Court, Wabasha, MN 55981

1-800-363-7846 | [www.medi-stim.com](http://www.medi-stim.com)