

CERVICAL

Part #POS1000

LUMBAR

Part #POS2000

POSTURE PUMP[®] DISC HYDRATOR[®]

Decompress, Shape & Lubricate™

Relieves Pain & Joint Stiffness
and Decreases Bulging Discs in the

← **Neck & Low Back** →

by Decompressing
Joints, Hydrating
Discs, and Aligning
the Vertebra into
their Naturally
Curved Shape.



Lightweight, Portable & Easy To Operate



Handheld pumps have automatic safety valves to limit inflation and allow for gentle release at proven therapeutic levels.



Deflated



Inflated



Deflated



Inflated

Cervical Disc Hydrator

- Shapes and gently stretches the cervical curve to restore range of motion and relieve neck pain.
- Forehead Strap: 3 adjustable snap settings and allows for removable chin section.
- Weight: 2.2 lbs
- Dimensions: 15" L x 14" W x 5" H
- Six Month Warranty

Elliptical Back Rocker

- Built-in Rocking Mechanism gently shapes, stretches and strengthens the lower back.
- 2 uniquely angled air cells create multi-vector force, hydrating discs and shaping the spine.
- Weight: 2.2 lbs
- Dimensions: 13" L x 8.5" W x 3.5" H
- Six Month Warranty



Composition (Both Models)

Inside: Exclusive bladder fabric 100% expandable "ballistic" nylon
Outside: Exclusive outside bladder stitching allows for greater expanding & flexibility.
100% neoprene and waterproof. **Frame:** Engineered high impact "ABS" made to last.



217 Industrial Ct., Wabasha, MN 55981
www.medi-stim.com | sales@medi-stim.com
1-800-363-7846